



Phadia[®] 250

Designed to aid in the diagnosis of allergic and autoimmune diseases

Test results with clinical relevance

Flexibility

- Automated processing of EliA[™] and ImmunoCAP[®] assays
- An allergen panel detecting sensitization to more than 650 allergens and including more than 100 allergen components
- An autoimmunity panel for evaluation of over 20 diseases
- All tests can be measured from one patient sample
- User-configurable reflex testing functionality

Capacity

- Throughput: 60 results/hour
- Continuous random access
- All reagents available onboard
- 28-day storage of calibration curves
- Onboard dilutions
- Connectable to Inpeco based Laboratory Automation Systems (LAS)
- Clustered set-up for increased capacity

Reliability

- Excellent consistency over time and between countries, systems, laboratories and users*
- Remote troubleshooting and support via Phadia LabCommunity
- Confirmation of accuracy and reliability with Quality Club

Ideal for medium-sized laboratories delivering 400–2000 results per week

*Based on Quality Club data

Phadia[®] 250

Proven performance in more than thousand laboratories all over the world!

- ❶ Stop solution and development solution, 2 bottles each.
- ❷ Loading tray for 38 EliA or ImmunoCAP carriers in addition to the positions in the storage chamber.
- ❸ Cooled storage chamber for 180 barcode-controlled EliA and ImmunoCAP carriers.
- ❹ Sample diluent (6 bottles) and dilution plates (2 × 96 wells) for onboard dilution of samples.
- ❺ Barcode reader, for identification of reagents.
- ❻ Cooled conjugate compartment (6 bottles).
- ❼ Strip tray (8 positions) or calibrator strips and/or curve control strips.
- ❽ LCD, logical and intuitive instrument management via touch screen.
- ❾ Integrated storage for wash and rinse solution, as well as solid and liquid waste containers.
- ❿ Connectable to fixed piping and drain, enabling automated filling of rinse solution and waste drainage.
- ⓫ Sample racks (5 positions), each holding up to 10 patient and quality control samples.

Not shown in picture:
Stand alone PC for Phadia Prime software, connectable to mainframe Laboratory Information Systems (LIS) and LAS.



ImmunoCAP assays available for Allergic diseases

- ImmunoCAP Specific IgE
- ImmunoCAP Total IgE
- ImmunoCAP Specific IgG
- ImmunoCAP Specific IgG4
- ImmunoCAP ECP
- ImmunoCAP Tryptase

EliA assays available for the following Autoimmune diseases

- Rheumatoid Arthritis
- Celiac Disease
- Inflammatory Bowel Diseases
- Vasculitis / Goodpasture Syndrome
- Connective Tissue Diseases
- Thyroid Diseases
- Autoimmune Thyroiditis
- Primary Biliary Cirrhosis

Automation for walk-away productivity



- Estimated hands-on time 10–15 minutes. Length of work shift is adaptable to your needs.
- First result available in only 100 minutes, then a result every 60 second.

Find out more at thermoscientific.com/phadia

© 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Apr. 21, 2017/01

Distributed in Canada by / Distribué au Canada par:
SOMAGEN[™]
DIAGNOSTICS
www.somagen.com 1.800.661.9993
*Some products may not be authorized for sale in Canada. Please contact Somagen for specific enquiries. Le vente de certains produits n'est pas autorisée au Canada. Si vous plaît veuillez contacter Somagen pour des demandes spécifiques.

ThermoFisher
SCIENTIFIC